AQUATOOLS	INSTRUCTION FOR USE FILTRAY point of use filter	EA2MR2 09/2025
EN FILTRAY point of use filter ES Filtros terminales FILTRAY PT Filtros terminais FILTRAY Manufacturer : AQUATOOLS – 11, rue Lobineau – 75006 PARIS France Main location site : AQUATOOLS – 18, rue du Maréchal Foch – 80130 FRIVILLE France US commercial subsidiary/Filial comercial EE. UU./Filial comercial EUA: AQUATOOLS Inc - 1 Broadway, 14th floor, Cambridge MA 02142 United States Tel: +1 516 464 1513		

1. Product compliance

Sterile* FILT'RAY Disposable Water Filters:

- Class II Medical Device in the US - Compliant with FDA 510(k) Clearance (510k Number: K191563)
Device are manufactured in compliance with the Quality System Regulation (21 CFR Part 820)
- Are **Class Is Medical Devices for Europe** that comply with European Regulation EU 2017/745 CE Marking year 2023 by notified body n° 0459.

**Sterile by Gamma irradiation in compliance with ISO 11137*

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- Are **Class I Medical Devices for Europe** that comply with European Regulation EU 2017/745 CE Marking year 2017.

Both comply with European Drinking Water Certifications (**ACS** - France, **KTW** - Germany, **WRAS** - United Kingdom).

Laboratory-based microbiological challenges of *Brevundimonas diminuta*, *Legionella pneumophila* and *Pseudomonas aeruginosa* on the FILT'RAY demonstrate a 7-log retention capacity.

2. References

Sterile* :

- FILT'RAY Shower head water filters: AT030711, AT030713, AT030714, AT030716.
- FILT'RAY Wall-mounted shower water filters: AT030777, AT030778, AT030779, AT030780.
- FILT'RAY Tap water filters with straight flow outlet: AT030717, AT030719, AT030747, AT030749.
- FILT'RAY Tap filters with shower flow outlet: AT030718, AT030720, AT030748, AT030750.
- FILT'RAY In-line water filters: AT030733, AT030734, AT030745, AT030746.

Non-sterile* :

- FILT'RAY Shower head water filters: AT030715, AT030723, AT030724, AT030725.
- FILT'RAY Wall-mounted shower water filters: AT030741, AT030742, AT030767, AT030768.
- FILT'RAY Tap water filters with straight flow outlet: AT030721, AT030727, AT030729, AT030730.
- FILT'RAY Tap water filters with shower flow outlet: AT030722, AT030728, AT030739, AT030740.
- FILT'RAY In-line water filters: AT030743, AT030744.

To guarantee the optimal performance of the FILT'RAY Disposable Water Filters, please ensure that you follow the instructions for installation (see) and read the warnings carefully.

**Sterile : Gamma sterilization during production process compliant to ISO 11137.*

3. Indication for use

- Point of use water filter - Disposable (single use)
- Available on version sterile (Gamma sterilization during production process compliant to ISO 11137) or non-sterile
- Filtration of potable water : 0.2 µm physical retention of microorganisms such as *Legionella Pneumophila* and *Pseudomonas Aeruginosa*.
- For up to 1, 2, 3 & 4 months
- They are intended to be installed in hospitals or medical care facilities (retirement homes, specialized centers), as well as in collective facilities such as leisure or sports centers, hotels.

4. Technical specifications

Conditions of use:

The performance of FILT'RAY filters is guaranteed when the pressure and temperature of the water in the network are respectively between 1 and 5 bar/41 °F and 140 °F. The measured outlet water flow rates are as follows:

Typical Flow Rate (gal/min)					
Product	Water Pressure (PSI)				
	14.5	29	43.5	58	72.5
FILT'RAY Shower Head Filter and Wall Mounted Shower Filter	1.6	2	2.4	2.7	3.1
FILT'RAY Faucet and In-Line Filter	1.1	1.3	1.6	1.9	2.1

** Flow rate measurements are determined with the filters installed on a faucet.*

Typical Flow Rate (L/min)					
Product	Water Pressure (Bar)				
	1	2	3	4	5
FILT'RAY Shower Head Filter and Wall Mounted Shower Filter	6	7.4	9	10.4	11.6
FILT'RAY Faucet and In-Line Filter	4	5	6	7.1	8

** Flow rate measurements are determined with the filters installed on a faucet.*

Exposure to chemical and thermal treatments:

The filter can tolerate the following treatments:

- Heat treatment of 70°C (158 °F) for 30 minutes (maximum cumulative).
- Chemical treatment: active chlorine.
- Filter 1M : 100 ppm for 1 hour.

- Filter 2M : 100 ppm for 2 hours.
- Filter 3M : 100 ppm for 3 hours.
- Filter 4M : 100 ppm for 4 hours.

Replace the filter if the chemical or thermal treatment exceeds these thresholds.

5. **Contraindication**

- Disposable Water Filters must not under any circumstances be used to produce water for **pharmaceutical purposes** (i.e., water that will be used in solutions intended to be administered orally, injected, or perfused).
- Disposable Water Filters are single-use (single use). They are sterilized by Gamma irradiation in compliance with ISO 11137. It does not mean that the filtered water is sterile. Do not re-sterilize or re-install on another point of use.

6. **Precautions for use**

Packaging and duration

- Check that the single packaging bag is not punctured or torn and that the filter is not damaged.
- Recommended Storage: Between 5 and 35°C (41 and 95°F).
- Check that the expiry date has not passed. If necessary, discard the filter because its safety and performance are no longer guaranteed. This creates a risk of cross-contamination.
- Avoid contact with strong solvents such as acetone, benzene, trichloroethylene, and perchloroethylene, or ammonia-based cleaners. Use of these disinfectants can damage microfilter membranes and result in discoloration and weakening of housing material.

Using

- Disposable Water Filters must not under any circumstances be used to produce water for **pharmaceutical purposes** (i.e., water that will be used in solutions intended to be administered orally, injected, or perfused).
- Disposable Water Filters are single-use (disposable). They cannot be sterilized or re-installed on another point of use.
- The filter does not change the potability of the water, and it does not modify its concentration of salts, metals, scale, chemicals, pesticides or dyes.
- Shower filters must always be repositioned on the support after use.

Clogging tolerance:

If the water in the drinking water network has an unusual rate of particles or microorganisms in suspension, it is possible to observe a decrease in the output flow rate before the theoretical decommissioning date (date of entry into service + theoretical duration of use).

When the flow rate becomes insufficient, replace the filter.

Warning: Clogging may occur before the theoretical discharge date.

Tip: In case of premature clogging, pre-filtration of the water system is recommended to reduce suspended solids and colloids.

7. **Installation instructions**

Preliminary preparation of the taps:

- The screw-on shower filter: it is screwed quickly clockwise and until the seal is closed on all shower hoses complying with standard NF EN 1113.
 - The length of the hose must be adapted to avoid any contact between the shower filter and the ground (70 cm/27.6 in.). The gasket provided must be installed to ensure sealing and prevent possible leaks. Each time the shower head is changed, it must be replaced.
 - The faucet filter: it must be fitted using the quick connectors recommended by AQUATOOLS. The quick connectors allow the filter to be connected to and disconnected from the faucet easily.
- Visit the website or contact us for a list of available fittings.

General Information:

- The faucet and wall-mounted shower filters can be connected to most plumbing fixtures using an adaptor (clip-on); contact us for a list of commercially available adaptors.

Installing:

1. Open the faucet to purge before installation to avoid early clogging.
2. Close or shut off the water supply to an electronic faucet.
3. When installing for the first time:
 - Faucet filter: place an adapter fitting.
 - Screw-on shower filter: remove the original shower head.
4. If replacing a worn filter: remove the filter.

5. Wash and disinfect your hands.

6. Partially open the unit packaging.
7. While keeping the filter in the bag to prevent contamination of the water outlet, screw the shower/hand spray filter to the end of the hose or clip the cartridge filter to the fitting adapter.
8. Remove the packaging.
9. Verify that the filter and fittings are tightly sealed by turning the water on for a few seconds after installation.
10. In case of leakage, repeat the procedure. Do not re-use a filter that has been connected.
11. Note the commissioning and the theoretical decommissioning date (commissioning date + theoretical duration of use) on the filter and on the traceability label.

Removal :

1. Close the water supply.
2. Release the used shower head filter by unscrewing it from the end of the hose or unclip the shower head filter or faucet filter from the adapter fitting.

3. Dispose of the filter as normal/household waste. Never re-use a filter that has been connected.

As the filter cannot be dismantled, the manufacturer cannot guarantee that the fibers are always in good condition.

Installation and use in a dirty or uncontrolled environment can be a source of retro-contamination.

Discard and replace any filter that appears visually soiled or contaminated.

However, use in a non-controlled environment can be a source of retro-contamination.

Any filter that is believed to be contaminated must be replaced, especially if the contaminant is visible.

Waste management:

Used FILT'RAY filters are classified as normal household waste. No special infectious risk.

The installation guide and the white carton can be recycled (paper/card).

Maintenance:

To track maintenance procedures, a removable adhesive label (located on the filtration device) can be used to record important information for your safety records (date of installation, date of removal, filter unique identification number, etc.).

This information is also included in the UDI (Unique Device Identification) which is a datamatrix with information such as (01)

GTIN: Global Trade Item Number, (10) Lot number, (17) expiration date, (21) Serial Number, (240) AQUATOOLS Part Number.

Manufacturing date, Expiration date and Lot number are also clearly identified on device labels and packaging.

The device label provides an additional space for the user to note clearly the date of first use and the date of disposal.

Duration of life is identified by a numerical symbol featuring the number of months (1, 2, 3, 4) of the use-life specification.

8. Manufacturing information

FILT'RAY filters are manufactured by:

Registered Office: **AQUATOOLS** | 11, rue Lobineau | 75006 PARIS | France

Main location site: **AQUATOOLS** | 18, rue du Maréchal Foch | 80130 FRIVILLE | France

EPA Est. No. 104636-FRA-1

Call us: +1 39 75 02 20

Distributed by US commercial subsidiary : AQUATOOLS Inc – 1 Broadway, 14th floor, Cambridge MA 02142 United States

Call us: +1 516 464 1513

Distributed by UK commercial subsidiary: AQUATOOLS HYGIENE LTD - Henderson House, Hithercroft Road, Wallingford, OX10

9DG UNITED KINGDOM

Call us: +1 491 521 282

1. For more information, visit our website: **INT:** www.aqua-tools.com

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

US: www.aquatools.us - **E-mail:** sales@aquatools.us

UK: www.aquatools.uk - **E-mail:** contact@aquatools.uk

2. The paper version of the notices remains available on request.

3. Any serious incident must be notified immediately to the manufacturer and to the competent authority of the Member State where the user and/or the patient resides.

9. Product label – Meaning of symbols

REF	LOT	
Catalog number	Batch Code	Expiration date
STERILE R		
Device sterilized by Gamma irradiation during production	Do not re-sterilize the device.	Do not use if the package is damaged
	GTIN	
Manufacturer's address	Global Trade Item Number	Consult the instructions for use
	MD	UDI
Sterile barrier system only	Medical Device	Unique Device Identifier
		
Do not reuse	Manufacturing date	Shower head filter
		
Tap filter	Shower flow	
SN		
Serial Number	Commissioning date	Decommissioning date
		Rx Only
Use-by date after commissioning 1M: 1 month - 2M: 2 months - 3M: 3 months - 4M: 4 months	Quantity	Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner
FILT'RAY	CE 0459	
Name of the range	Meets general safety and performance requirements	In-line filter
		
Flow direction of in-line filters		

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