


# FILT'RAY Scan

INFECTION CONTROL  
SOLUTIONS

## Secure Point-of-Use Water Filtration Management

**TRACK FILT'RAY<sup>2G</sup> MAINTENANCE ACTIONS**

Android and iOS apps





1

View point-of-use water filters

2

Uninterrupted access to functionalities

3

Datamatrix scan to add or replace filters

**AQUATOOLS** mobile app oversees and controls the maintenance of **FILT'RAY<sup>26</sup>** showerhead and faucet Point-of-Use Filters in facilities and locations open to the public.

Sound hygiene and safety in healthcare establishments requires strict Infection Risks Management at water points-of-use combined with rigorously disciplined tracking of filter maintenance and replacement.

Depending on type, **FILT'RAY<sup>26</sup> filters** are designed for optimal safety when used from **1 - 2 months** or **up to 3 months**. Filters should then be replaced.

**FILT'RAY Scan** is an app securely accessible online by tablet or smartphone.

First, Water Point-of-Use is entered and identified (by room number or other). Then, **FILT'RAY Scan** captures and displays the number of days left before filter replacement. It also sets up and monitors inspection schedules, and records software user access privileges and the history of maintenance actions.

**FILT'RAY Scan's** Maintenance Monitoring capability makes it ideal in every healthcare facility for recording and planning maintenance of point-of-use water filtration.

## BENEFITS

Secure personalisable platform  
E-mail or SMS messaging to prompt timely filter replacement

Comprehensive traceability of **FILT'RAY<sup>26</sup>** filter replacement

Scanned capture and tracking of filter identification data

Show number of days left until filter replacement

Schedule and track completion of maintenance actions

Determine end-user privileges

Draft and dispatch reports in Excel and PDF formats

